



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Fractional CO2 Vaginal Laser (New)

- New,
- Made in South Korea,
- Advantages:
 - Maximum power: 35 W,
 - Individual adjustment of power and pulse,
 - 4 surgical handpieces, fractional scanner, endovaginal scanner,
 - Mass storage capacity,
 - Max power: 1050 mJ at 3x overlapping,
- Indication:
 - Vaginal contraction,
 - Improvement of urinary incontinence,
 - Improvement in vaginal dryness,
 - Improved vaginitis,
 - Improves sexual satisfaction,
- Technical Specifications:
 - Laser Type: Laser CO2 Sealed Off, 10600nm,
 - Touch screen,
 - Ultra Pulse Mode:
 - Surgical:
 - Continuous Wave Mode: 1 W ~ 35 W,
 - Ultra Puls Mode: 0,5 - 35 W, Pulse width: 90 - 900 us,
 - Super Puls Mode: 0,5- 35 W, Pulse width: 10 ms - 20 ms,
 - Basic Puls Mode: 0,5 - 35 W, Pulse width: 3 - 100 ms,
 - Fractional:
 - Power: 0,2 - 350 mJ,
 - Pulse width: 90 us - 2 ms,
 - Overlapping: 1-5x,
 - Endovaginal:
 - Power: 0,2 - 1000 mJ, beauty impuls,
 - Pulse width: 90 us - 1000 mJ,
 - Overlapping: 1-3x,
 - Aiming Beam: 6 mW, 650 nm red diode laser,
 - Delivery system: 7-joint, fixed mirror articulated arm,
 - Spot size: ≥ 120 microns,
 - Scanning area: 20x20 mm
 - Scanning area of endovaginal handpiece: 10,8 x 10,8 mm
 - Scanning direction: sequence, diverge, random,
 - Scanning shape: square, rectangle, circle, triangle, hexagon,
 - Dot pitch (density): 0,2- 2,5mm each 0,1 mm,
 - Dimensions: 368 x 377 x 1050,
 - Power input: 100 ~ 230 VA, 50/60 Hz,
 - Maximum power consumption: 1500 VA,
- Technical Passport,
- 24-month warranty

