



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Fractional CO2 Laser S-CO2 Galvano (New)

- New,
- Made in South Korea,
- Used in major specializations:
 - Aesthetic Medicine,
 - Dermatology,
 - Surgery,
 - Laryngology,
 - Proctology,
 - Gynaecology,
- Technical Specifications:
 - Wavelength - 10600 nm,
 - Beam Mode Quality - TEM00,
 - Laser power output - 1-35 W,
 - Spot size - 120-300 μ m,
 - Aiming Beam - 3,5 mW, 650 nm,
 - Cooling System - liquid,
 - Ultra Pulse Mode:
 - Pulse width: 100 - 2000 Us,
 - Repeat delay: single, 0,1 - 500 ms,
 - Power: 1 - 35 W,
 - Super Pulse Mode:
 - Pulse width: 1000 - 2000 Us,
 - Repeat delay: single, 0,1 - 500 ms,
 - Power: 1 - 35 W,
 - Basic Puls Mode:
 - Pulse width: 3 - 100 ms,
 - Repeat delay: single, 0,1 - 500 ms,
 - Power: 1 - 35 W,
 - Continuous Wave Mode: 1 - 35 W,
 - Maximum power consumption- 1500 VA,
 - Power input: 240 V, 50/60 Hz,
 - Dimensions: 360 x 430 x 1200
 - Weight: 54 kg,
- Scanner technical specifications:
 - Power: 1-35 W,
 - Scanning area: from 0,4 x 0,4 to 20 x 20 mm,
 - Scanning shape: square, rectangle, circle, triangle, hexagon,
 - Scanning direction: sequence, diverge, random,
 - Dot pitch (density): 0,1 - 2,5 mm,
- Technical Passport,
- 24-month warranty

