



**PRAISTON Sp. z o.o.**  
 ul. Górowska 32  
 64-100 Leszno  
 Polska

Biuro: +48 65 527 01 67  
 Serwis: +48 782 844 000  
 Fax: +48 65 527 01 67  
 E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
 Krzysztof Wybieralski  
 +48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



**N.I.L. EVA sp6 (New)**

EVA sp®6 is composed by a control box (with an entirely digitalized display) regulating the compressed air supply to the Lipomatic® which allows the Nutational Infrasonic Liposculpture®, a high-speed aspirator and a quantitative infiltration system. EVA sp®6 is a complete system. The quantitative infiltration system allows to control precisely the infiltrated quantity in different predefined zones. The system of pistons ensures an homogeneous infiltration at a constant pressure, whatever the type of fat, even fibrous.

System Composition:

- 1x Lipomatic® 3,
- 2x Pedals,
- 5x Sewer jetter,
- 5x Air tubings,
- 5x Perfusion lines,
- 5x Suction liner Medivac Flex 3L,
- 5x Infiltration tubings,
- BAXTERS (Set):
- 1x Power Supply Unit,
- 1x Ragulator Key,
- 1x Cleaning Set,
- 1x Lipomatic® Oil (125 ml),
  
- 2x Canister jars,
- Complete Set of Cannulas:
  - 1x Mercedes 4x 300,
  - 1x 3G Accelerator 3x 250,
  - 1x 3G Accelerator 3x 150,
  - 1x Rebelo 3,5x 250,
  - 1x Infiltration 2.5x 200,
  - 1x Infiltration 2.5x 300.

If you have any questions, please do not hesitate to contact us!  
 In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

- WE SPEAK ENGLISH
- HABLAMOS ESPAÑOL
- WIR SPRECHEN DEUTSCH

