



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **UZUMCU OM-6N Operating OR table (New)**

- New,
- Manufactured in Turkey,
- Features and advantages of OM-6N table:
  - The table is made of stainless steel (AISI 316L),
  - The mattress is covered with a special antimicrobial coating that protects the surface from bacterial growth while not allowing blood and other fluids present in the operating room to penetrate,
  - The X-ray-permeable top plate provides high-quality images in all parts of the table,
  - The table base has a T-shaped design, in addition, with a 300 mm sliding feature, the table provides easy access and approach to the table for instruments such as a C-arm X-ray renderer and operating microscope,
  - Accessory rails around the operating table make it easy to position all the equipment to be used during surgery and create a comfortable working environment,
- Electrically adjustable:
  - Height,
  - Trendelenburg position,
  - Anti-Trendelenburg Position,
  - Backrest,
- Mechanical adjustment:
  - Headrest,
  - Leg sections: up/down,
  - Leg sections: left/right,
- Technical specifications:
  - Length x width : 2200 x 610mm,
  - Height (min) : 800mm,
  - Height (max) : 1200mm,
  - Backrest adjustment : +90° / -40°,
  - Trendelenburg position : +25°,
  - Anti-Trendelenburg position : -25°,
  - Lateral tilt : +20° / -20°,
  - Leg rest adjustment (up/down) : +65° / -80°,
  - adjustment of leg rest inversion : 90°,
  - Headrest section adjustment : +30° / -80°,
  - Flex/Reflex position : +220° / 90°,
  - Upholstery thickness : 50 mm,
  - Number of segments : 6,
  - Wheel diameter: front wheel Ø80x30mm,
  - Double rear wheel: Ø125x32mm,
  - Number of wheels: 3,
  - Maximum patient weight: 300 kg,
  - Position for kidney surgery: +250° / 12 cm,
  - Table base construction: T-type,
  - Control controller: handheld remote control,
  - Memory function: ability to store 5 different positions,
  - One-touch position reset: YES,
- Possibility to retrofit the table with accessories visible in the last 2 photos (for an extra charge),
- Valid technical passport issued,
- Warranty 24 months,
- Possibility of financing (Poland only): Installments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPAÑOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANÇAIS

WIJ SPREKEN DUTCH

Мы говорим по русски

Мы говорим по украински

