



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **DERMATIC 1 - Meso Gun (New)**

Dermatic®1 is a high technology gun for mesotherapy functioning with compressed air. The Dermatic® principle consists in using compressed air that provides a high-injection speed. This ensures a painless treatment for the patient and a very high precision for the practitioner.

### **Advantages for the surgeon:**

- No loss of product,
- Excellent ergonomy,
- Regulation of the injected quantity,
- Constant safety,
- Easy regulation of the parameters,
- Manual or automatic mode,

### **Advantages for the patient:**

- Less painful for the patient,
- Speed of burst,
- Injection depth,
- Quantity of liquid to inject,
- Constant safety,

### **System Composition:**

- 1x Dermatic®1 (Meso Gun),
- 1x Suitcase,
- 20x Tips + one time catheters (with syringes),
- 1x Tip (sterile),
- 20x Syringes (10 ml),
- 1x Syringe Adapter 10 ml.

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPAÑOL

WIR SPRECHEN DEUTSCH

