



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Medical Display BEACON C24S+ (New)

- New,
- 2MP Color Medical Display,
- Features:
  - MULTI,
  - BSS,
  - ALAS,
  - DICOM,
  - COLOR,
- Application:
  - CR/DR,
  - CT,
  - MRI,
  - PET / CT,
  - X-RAY,
  - DSA,
  - DSI,
  - 3D,
  - Nuclear Medicine,
  - CAD medical software,
  - Image reconstruction,
- Technical Specifications:
  - panel Size: 21,3",
  - Support Color: 1,07 Billions,
  - Resolution: 1600 x 1200,
  - Pixel Pitch: 0,270 x 0,270 mm,
  - Brightness: 440 cd / m2,
  - Contrast: 1500 : 1,
  - Angle of View: 178°,
  - Response Time: 20 ms,
  - Input:
    - DVI-D: 24 pin,
    - DisplayPort: 20 pin,
    - VGA: D-Sub 15 pin,
  - Dimensions: 369 (W) x 507 ÷ 627 (H) x 235 (D) mm,
  - Weight: 9,2 kg,
  - Standard Mounting: VESA 100 x 100 mm,
  - Power supply: DC 12 V,
  - Max power consumption: 60 W,
  - Power Saving: EPA, VESA,
  - MTBF: > 50 000h,
  - Certificates: EN 60601, CCC, FCC, CE, UL,
- Technical Passport,
- 36-month warranty

