



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **Medical Display BEACON C14S (New)**

- New,
- 1MP Color Medical Display,
- Features:
  - MULTI,
  - PERSPEKTIVE,
  - BSS,
  - DICOM,
- Application:
  - DSA,
  - DSI,
  - CT,
  - MRI,
- Technical Specifications:
  - Panel Size: 19",
  - Support Color: 16,7 Billions,
  - Resolution: 1280 x 1024,
  - Pixel Pitch: 0,294 x 0,294 mm,
  - Brightness: 300 cd/m<sup>2</sup>,
  - Contrast: 1000 : 1,
  - Angle of View: 178°,
  - Response Time: 14 ms,
  - Input:
    - DVI-D: 24 pin,
    - DisplayPort: 20 pin,
    - VGA: D-Sub 15 pin,
  - Dimensions: 425 (W) x 449 ÷ 523 (H) x 225 (D) mm,
  - Weight: 6,3 kg,
  - Standard Mounting: VESA 100 x 100 mm,
  - Power supply: AC 230 V,
  - Max power consumption: 50 W,
  - Power saving: EPA, VESA,
  - MTBF: > 50 000 h
  - Certificates: CCC, CB, FCC, CE, cTUVus,
- Technical Passport,
- 36-month warranty

