



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **Medical Display BEACON C22WT (New)**

- New,
- 2MP Color Medical Display,
- Features:
  - SUPPORT TOUCH,
  - PROTECTIVE,
  - DISPLAY,
  - BSS,
  - DICOM,
- Application:
  - Outpatient Service,
  - ICU (Intensive Care Unit),
  - Clinic Department,
  - Community Medical Service,
- Technical Specifications:
  - Panel Size: 21,5",
  - Support Color: 16,7 Billions
  - Resolution: 1920 x 1080,
  - Pixel Pitch: 0,248 x 0,248 mm,
  - Brightness: 250 cd/m<sup>2</sup>,
  - Contrast: 1000 : 1,
  - Angle of View: 178°,
  - Response Time: 14 ms,
  - Input:
    - DVI-D: 24 pin,
    - DisplayPort: 20 pin,
    - VGA: D-Sub 15 pin,
    - Audio 3,5 mm stereo jack,
  - Speakers: 5 W\*2,
  - Dimensions: 523 (W) x 376 (H) x 179 (D) mm,
  - Weight: 6 kg,
  - Standard Mounting: VESA 100 x 100 mm,
  - Power supply : AC 100 - 240 V / DC 12 V,
  - Max power consumption: 40 W,
  - Power Saving: EPA, VESA,
  - MTBF: > 50 000h.
  - Certificates: EN 60601, CCC, CB, FCC, CE, cTUVus,
- Technical Passport,
- 36-month warranty

