



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Display BEACON S421 (New)

- New,
- FHD Surgical Color LCD Display,
- Features:
 - MULTI,
 - RESOLUTION FULL HD,
 - DUAL,
 - INPUT/OUTPUT,
 - LED,
 - COLOR,
 - BSS,
- Application:
 - Endoscopy,
 - Operating room,
 - Gastroscopy/Laparoscopy,
 - ICU (Intensive Care Unit),
- Technical Specifications:
 - Panel Size: 42",
 - Support Color: 1,07 Billions,
 - Resolution: 1920 x 1080,
 - Pixel Pitch: 0,4833 x 0,4833 mm,
 - Brightness: 700 cd/m²,
 - Contrast: 1300 : 1,
 - Angle of View: 178°,
 - Response Time: 8 ms,
 - Input:
 - DVI-D: 24 pin,
 - VGA: D-Sub 15 pin,
 - CVBS: BNC,
 - S-VIDEO: DIN4,
 - YPbPr: BNC*3,
 - RGBS: BNC*4,
 - HD-SDI,
 - 3G-SDI ,
 - FC,
 - Output:
 - DVI-D,
 - SDI,
 - DC 5V,
 - Dimensions: 1003 (W) x 597 (H) x 65 (D) mm,
 - Weight: 45 kg,
 - Standard Mounting: VESA 400 x 400 mm,
 - Power supply: AC 230 V,
 - Max power supply 120 W,
 - Power saving: EPA, VESA,
 - MTBF: > 60 000 h,
 - Certificates: EN 60601, CCC, FCC, CE, UL,
- Technical Passport,
- 36-month warranty

