



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Autoclave STATIM 5000 S (no registration) (New)**

- New,
- Made in Canada,
- Features:
  - Ultra fast 8:45 minute unwrapped cycle,
  - 15:30 minute wrapped cycle,
  - Exceptional instrument turnover,
  - 15 minute Dri-Tec drying technology,
  - Single use water,
  - Aseptic transfer of instrument,
  - Reduced instrument investment,
  - Gentle on instruments,
  - EN13060 compliant,
  - Printer - optional,
- Technical Data:
  - Class S,
  - Power input: 208 - 240 V,
  - Power: 1300 W,
  - Decibel rating - 56 dB,
  - Cycle registration- optional Data Logger,
  - Hollow wrapped 134 C - 18 mins,
  - Hollow unwrapped 134 C - 9 mins,
  - Rubber & Plastic 121 C - 24 mins,
  - Drying cycle minimum 60 mins,
    - Wrapped 82 mins,
    - Unwrapped - 45 mins,
    - A + B,
  - Cassette dimensions - 37 cm / 18 cm / 6,5 cm,
  - Cassette capacity - 5,1 l,
  - Reservoir capacity: 4 l,
  - Dimensions: 55 cm x 41 cm x 19 cm,
  - Menu in Polish,
  - Weight - 33 kg,
- Technical Passport,
- 24-month warranty

