



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Autoclave STATIM 5000 G4 WIFI (Data Logger built-in) (New)**

New,  
Canadian production,  
Power supply: 230 V,  
Frequency: 50/60 Hz,  
Power: 1300 W,  
Selectable menus in: Polish, English, Czech, Danish, Spanish, French, Hungarian, Italian, Korean, Dutch, Portuguese, Romanian, Russian, Ukrainian,  
Advantages:  
Colour touchscreen LCD display,  
Advanced G4 operating system,  
Easy to use,  
Built-in Data Logger,  
USB flash drive recording,  
LAN network interface,  
Energy saving,  
Fast sterilization cycles - up to 9 min. for unpacked instruments,  
Low consumption of distilled water,  
Simple construction facilitating service and periodical maintenance,  
Compliance with EN13060 for sterilisation of type A and B loads,  
Technical data:  
Class S,  
Operating noise - 57 dB,  
Cycle recording - built-in Data Logger,  
Cycle for packed 134°C - 18 min,  
Cycle time for unpacked 134°C - 9 min,  
Plastics and rubber 121°C - 24 min,  
Drying - minimum 60 min,  
Sterilisation and drying of packaged - approx. 82 min,  
Sterilisation and drying of non-packaged - approx. 45 min,  
Sterilisation of batches of type - A and B,  
Dimensions of the cassette (L/W/D) - 38 x 18 x 8 cm,  
Capacity of the cassette - 5.1 l,  
Water tank capacity - 4 l,  
LAN operation : WiFi or network cable RJ45,  
Equipped with RS232 port,  
Pre-vacuum: thermodynamic,  
Drying: thermal,  
Included:  
STATIM 5000-S cassette with shelf with drying accelerator plates,  
STAT-DRI Plus 59ml with atomiser,  
Test sample - Process Challenge Device,  
Sterilisation emulsifiers - 20 pieces,  
Waste bottle,  
Outlet pipe,  
Fittings for tube attachment,  
USB memory stick 16GB,  
USB memory stick with user manual in: Polish, Czech, Danish, Spanish, French, Hungarian, Italian, Korean, Dutch, Portuguese, Romanian, Russian, Ukrainian,  
User manual in English,  
Power cord,  
Dimensions: 58 x 41 x 20 cm,  
Set weight: 40 kg,  
Holds a current Technical Passport,  
Warranty: 12 months,  
Finance options (Poland only): Instalments, Leasing, Loan

