



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Autoclave STATIM 5000 G4 ENDO (Data Logger) (New)

- New,
- Made in Canada,
- Features:
 - Monitoring real-time cycle data,
 - Instant data logging,
 - Proactive maintenance management,
 - Automatic system updates to your G4 machines,
 - Single use water,
 - Aseptic transfer of instrument,
 - Reduced instrument investment,
 - Gentle on instruments,
 - EN13060 compliant,
- Technical Data:
 - Class S,
 - Power input: 208 - 240 V,
 - Power: 1300 W,
 - Decibel rating - 56 dB,
 - Cycle registration- Data Logger,
 - Hollow wrapped 134 C - 18 mins,
 - Hollow unwrapped 134 C - 9 mins,
 - Rubber & Plastic 121 C - 24 mins,
 - Drying cycle minimum 60 mins,
 - Wrapped 82 mins,
 - Unwrapped - 45 mins,
 - A + B,
 - Cassette dimensions - 37 cm / 18 cm / 6,5 cm,
 - Cassette capacity - 5,1 l,
 - Reservoir capacity: 4 l,
 - Dimensions: 55 cm x 41 cm x 19 cm,
 - Menu in English,
 - Weight - 33 kg,
- Technical Passport,
- 24-month warranty

