



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Sealing Device EUROSEAL 2001 Plus (New)

- New,
- Made in Italy,
- Features:
  - The sealing lever at the front makes packing easier,
  - The roll-locking bar keeps the roll stationary inside the sealing area at all times,
  - It is extremely ease to use and it limits the possibility of mistakes,
  - The high-quality materials and solid design make it indestructible,
  - The compact size allows it to be placed on a work surface or flush-mounted,
  - The smooth surfaces allow quick and thorough cleaning,
- Power supply - 230 V,
- Power - 110 W,
- Sleevs handle - 39 cm,
- Sealing width - 12 mm,
- Dimensions: 48 cm x 39 cm x 19 cm,
- Weight - 7,4 kg,
- Technical Passport,
- 12-month warranty

