



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Endoscopy table AION R02 (New)

New,  
Made in Poland,  
It has a simple and stable frame made of powder-coated steel profiles,  
Adjustable height and three-section sleep surface with standard adjustable headrest with cut-out for face and plug,  
Step height adjustment of the table,  
Headrest is adjustable with a gas spring, whereas the footrest is adjustable with a rastomat mechanism,  
It is possible to order upholstery in a colour of your choice free of charge, except for "PRESTIGE" colours which are available for an extra charge,

Additional options are available (at extra charge):

Underlay hanger,

Bottom shelf,

Armrests,

Castors with brake,

Technical parameters:

Length 190 cm,

Width 55/70 cm,

Height 55/70 cm,

empty weight 37 kg,

Maximum load 150 kg,

Headrest angle adjustment -700 to +450,

Adjustment of footrest angle 00 to 300,

It has a valid Technical Passport

Warranty on structure 24 months,

Warranty on hydraulic accessories 12 months,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

