



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Pace Tech SpO2 Intermediary Cable (New)

- New,
- Made in Uruguay,
- Sensor connection / device connection:
 - DB9 / REDEL 5 PIN (BIOOX.PT.IC2),
 - DB9 / REDEL 6 PIN (BIOOX.PTR.IC2),
- **Cable Length 2,0 m** (cable extension available on demand),
- Compatible with:
 - Generra: 710, 200, 300, 400, 810, 530; Pace Tech: 520,
 - Vitalmax 4000, Vitalmax 4100 (1st version),
- 12 months warranty,

