



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## Temperature Probe - Air Shields - Hill Rom (New)

- New,
- Made in Uruguay,
- Options:
  - Adult skin,
  - Pediatric skin,
  - Air,
- **Cable Length 1,5 m** (cable extension available on demand),
- Compatible with:
  - C-100, C-100QT, C-200, C-450, C-450QT, 550QT; Warmer Systems: 7810, 7820, 7830, 7840,
  - Isolette C2000, Versalet 7700,
  - C-450,
  - C-100, C-100QT, C-200, C-450,
- 12 months warranty,

