



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Temperature Probe - Dräger-Siemens (New)

- New,
- Made in Uruguay,
- Options:
 - **Dräger 8000 Series**(Cable length 1,5 m):
 - Adult Skin (BIOTEMP.2001PT),
 - Pediatric Skin (BIOTEMP.2001PTP),
 - **Dräger Caleo Series**(Cable length 1,5 m):
 - Adult Skin (BIOTEMP.2200PT),
 - Pediatric Skin (BIOTEMP.2200PTP),
 - **Infinity Vista, Delta, Gamma, Kappa**(Cable length 3,0 m):
 - Adult skin (BIOTEMP.1100PT),
 - Pediatric skin (BIOTEMP.1100PTP),
 - Tympanic (BIOTEMP.1100T),
 - Adult esophageal/rectal (BIOTEMP.1100RA),
 - Pediatric esophageal/rectal (BIOTEMP.1100RP),
- **Cable extension available on demand,**
- 12 months warranty,

