



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Temperature Probe - GE Datex (New)

- New,
- Made in Uruguay,
- Options:
 - **All Ohmeda Care Plus devices** (Cable length 1,5 m):
 - Air (BIOTEMP.300A3),
 - **All Ohmeda Giraffe devices** (Cable length 1,5 m):
 - Adult skin (BIOTEMP.1700PT),
 - **All GE-Datex devices** (Cable length 3,0 m):
 - Adult skin (BIOTEMP.400PT),
 - Pediatric skin (BIOTEMP.400PTP),
 - Tympanic (BIOTEMP.400T),
 - Adult esophageal/rectal (BIOTEMP.400RA),
 - Pediatric esophageal/rectal (BIOTEMP.400RP),
- **Cable extension available on demand,**
- 12 months warranty

