



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67
 E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006
kw@praiston.pl



Temperature Probe - Mindray (New)

- New,
- Made in Uruguay,
- Options:
 - **T5 / T8 series, iPM-9800** (Cable Length 3,0 m):
 - Adult Skin (BIOTEMP.1400PT),
 - Pediatric Skin (BIOTEMP.1400PTP),
 - Timpanic (BIOTEMP.1400T),
 - Adult esophageal/rectal (BIOTEMP.1400RA),
 - Pediatric esophageal/rectal (BIOTEMP.1400RP),
 - **All Mindray devices with ¼ plug in, (Cable Length 3,0 m):**
 - Adult skin (BIOTEMP.400PT),
 - Pediatric skin (BIOTEMP.400PTP),
 - Timpanic (BIOTEMP.400T),
 - Adult esophageal/rectal (BIOTEMP.400RA),
 - Pediatric esophageal/rectal (BIOTEMP.400RP),
- **Cable extension available on demand,**
- 12 months warranty,

