



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Temperature Probe - Philips HP (New)

- New,
- Made in Uruguay,
- Options:
  - Adult skin (BIOTEMP.1200PT),
  - Pediatric skin (BIOTEMP.1200PTP),
  - Tympanic (BIOTEMP.1200T),
  - Adult esophageal/rectal (BIOTEMP.1200RA),
  - Pediatric esophageal/rectal (BIOTEMP.1200RP),
- **Cable Length 3,0 m** (cable extension available on demand),
- Compatible with:
  - M1029
- 12 months warranty,



Model	Accessories	Price
BIOTEMP.1200PT	BIOTEMP.1200PT	12000000
BIOTEMP.1200PTP	BIOTEMP.1200PTP	12000000
BIOTEMP.1200T	BIOTEMP.1200T	12000000
BIOTEMP.1200RA	BIOTEMP.1200RA	12000000
BIOTEMP.1200RP	BIOTEMP.1200RP	12000000