



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Temperature Probe - YSI 700 (New)

- New,
- Made in Uruguay,
- Options:
  - Adult Skin (BIOTEMP.700PT),
  - Adult esophageal/rectal (BIOTEMP.700RA),
  - Pediatric esophageal/rectal (BIOTEMP.700RP),
- **Cable Length 3,0 m** (cable extension available on demand),
- Compatible with:
  - All YSI 700 devices,
- 12 months warranty,

