



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **Operating table Uzumcu OPERA OP-2P (New)**

New,  
Turkish production,  
Electrically adjustable:  
height,  
Trendelenburg position,  
anti-Trendelenburg position,  
backrest,  
side tilt,  
translation,  
Mechanical adjustment:  
headrest,  
leg sections: up/down  
leg sections: left/right,  
central locking system,  
Technical specifications:  
the modular design is applicable for various surgical operations with high stability,  
thanks to the wide range of additional options the table is applicable for all kinds of surgical operations,  
4 large castors (125 mm) made of ABS plastic ensure mobility and easy manoeuvring of the table,  
table made of stainless steel (AISI 316L),  
table top is made of material transparent for X-rays,  
The mattress is covered with a special antibacterial coating that protects the surface against bacterial growth, at the same time preventing the penetration of blood and other fluids present in the operating room,  
length x width: 2130 x 610 mm  
height (min): 810 mm,  
height (max) : 1210 mm,  
travel: 300 mm  
backrest adjustment : +90° / -40°,  
Trendelenburg position : +25°,  
anti-Trendelenburg position : -25°,  
lateral tilting : +20° / -20°,  
adjustment of backrest for legs (up/down): +60° / -95°,  
adjustment of backrest angle: 80°,  
maximum load capacity: 350 kg (flat position) 200 kg (tilting position),  
weight of the table 200 kg,  
Possibility to add accessories shown in the last 3 photos (extra charge),  
It has a valid Technical Passport,  
Warranty: 24 months,  
Possibility of financing (only Poland): instalments, leasing, loan



PRAYSTON