



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **Gynaecological Patient Bed TEYCO MED ELLY04 (New)**

- New,
- Made in Italy,
- Equipped with 3 powerful actuators,
- Made entirely of steel painted with epoxy,
- Electrical Movements:
  - Height adjustment,
  - Trendelenburg,
  - Back Rest,
  - Leg Plates,
- Stainless Steel Irrigation bowl,
- Technical Specifications:
  - Power input: 220/230 V - 50 Hz,
  - Power: 24 V,
  - Height: 6000 N,
  - Back Rest: 3000 N,
  - Maximum Load: 180 kg,
  - Maximum Power: 120 VA,
  - Fuse: 1,6 + 6 A,
  - IP44 protection standard,
  - Weight: 90 kg,
  - Base Length: 800 mm,
  - Sitting height max: 840 mm,
  - Sitting height min: 480 mm,
- **Additional Equipment (extra paid):**
  - Telescopic pillar (additional or exchanged for leg plates),
  - Colposcopy support,
  - Headrest,
  - Paper dispenser,
  - Electrical Back Rest,
  - Electrical Leg Plates,
  - IV Pole,
  - Knee Clutch,
  - Casters with break system,
  - Knee Clutch cover,
  - Metal telescopic pillar model RING,
  - Chrome Base.
  - Transparent Glass Base: 19 mm.
- Technical Passport,
- 24-month warranty

