



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67
E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006
kw@praiston.pl



Examination Table TEYCO MED EXPEDIA E01 (Expedia top) (New)

- New,
- Made in Italy,
- Equipped with 2 powerful actuators,
- Folded footstool,
- Sitting: adjustable angle 45 ° left/right,
- Electrical Movements: Height, Back rest and footstool,
- Technical Specifications:
 - Power Input: 230 V - 50 Hz,,
 - Power: 24 V,
 - Height: 6000 N,
 - Back rest/Footstool: 6000 N,
 - Max Load: 180 kg,
 - Max power: 120 VA,
 - Fuse: 1,6 + 6 A,
 - IP44 protection standard,
 - Weight: 75 kg,
 - Base length: 580 x 740 mm,
 - Sitting height max: 980 mm,
 - Sitting height min: 650 mm,
- **Additional Equipment (extra paid):**
 - Headrest,
 - IV Pole,
 - Casters with breaks, metal frame,
- Technical Passport,
- 24-month warranty

