



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67  
E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **Examination Table TEYCO MED EXPEDIA E01 (Expedia top) (New)**

- New,
- Made in Italy,
- Equipped with 2 powerful actuators,
- Folded footstool,
- Sitting: adjustable angle 45 ° left/right,
- Electrical Movements: Height, Back rest and footstool,
- Technical Specifications:
  - Power Input: 230 V - 50 Hz,,
  - Power: 24 V,
  - Height: 6000 N,
  - Back rest/Footstool: 6000 N,
  - Max Load: 180 kg,
  - Max power: 120 VA,
  - Fuse: 1,6 + 6 A,
  - IP44 protection standard,
  - Weight: 75 kg,
  - Base length: 580 x 740 mm,
  - Sitting height max: 980 mm,
  - Sitting height min: 650 mm,
- **Additional Equipment (extra paid):**
  - Headrest,
  - IV Pole,
  - Casters with breaks, metal frame,
- Technical Passport,
- 24-month warranty

