



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67
 E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006
kw@praiston.pl



Phlebotomy Chair TEYCO MED LORY Y01 (New)

- New
- Made in Italy
- Equipped with 2 powerful actuators,
- Electrical Movements:
 - Height adjustment,
 - Back Rest,
- Technical Data:
 - Power supply: 230V - 50 Hz,
 - Power: 24 V,
 - Back Rest/Foot Rest: 6000 N,
 - Gas Spring: 800 N,
 - Max Load: 180 kg,
 - Max Power: 120 VA,
 - Fuse: 1,6 + 6 A,
 - IP44 protection standard,
 - Weight: 75 kg,
 - Base Length: 580 x 740 mm,
 - Sitting height max: 980 mm,
 - Sitting height min: 650 mm,
- **Additional Equipment (paid):**
 - Headrest,
 - Paper dispenser,
 - IV Pole,
 - 4 gummy wheels with metal frame + breaks,
- Technical Passport,
- 24 months warranty

