



**PRAISTON Sp. z o.o.**  
ul. Górska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Phlebotomy Chair TEYCO MED LORY Y01 (New)

- New
- Made in Italy
- Equipped with 2 powerful actuators,
- Electrical Movements:
  - Height adjustment,
  - Back Rest,
- Technical Data:
  - Power supply: 230V - 50 Hz,
  - Power: 24 V,
  - Back Rest/Foot Rest: 6000 N,
  - Gas Spring: 800 N,
  - Max Load: 180 kg,
  - Max Power: 120 VA,
  - Fuse: 1,6 + 6 A,
  - IP44 protection standard,
  - Weight: 75 kg,
  - Base Length: 580 x 740 mm,
  - Sitting height max: 980 mm,
  - Sitting height min: 650 mm,
- **Additional Equipment (paid):**
  - Headrest,
  - Paper dispenser,
  - IV Pole,
  - 4 gummy wheels with metal frame + breaks,
- Technical Passport,
- 24 months warranty

