



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Chair TEYCO MED RENE R01 (New)

- New,
- Made in Italy,
- Equipped with 1 powerful actuators,
- Electrical Movements:
 - Height adjustment,
- Height/Foot Rest/Trendelenburg regulated by gas spring,
- Technical Specifications:
 - Power input: 230 V – 50 Hz,
 - Power: 24 V,
 - Height: 6000 N,
 - Maximum load: 180 kg,
 - Maximum Power: 120 VA,
 - Fuse: 1,6 + 6 A,
 - IP44 protection standard,
 - Weight: 90 kg,
 - Base length: 800 mm,
 - Maximum height of the foot rest: 950 – 1000 mm,
 - Sitting height max: 950 mm,
 - Sitting height min: 500 mm,
- **Additional Equipment (extra paid):**
 - Headrest,
 - Paper dispenser,
 - Casters - powder coated metal,
 - Chrome base,
 - IV Pole,
 - Enamelled base,
 - Transparent crystal glass,
 - Base from transparent crystal glass: 19 mm,
- Technical Passport,
- 24-month warranty

