



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Examination Chair TEYCO MED RENE R03 (New)**

- New,
- Made in Italy,
- Equipped with 3 powerful actuators,
- Electrical Movements:
  - Height adjustment,
  - Back Rest,
  - Trendelenburg,
- Foot Rest regulated by gas spring,
- Technical Specifications:
  - Power input: 230 V - 50 Hz,
  - Power: 24 V,
  - Height: 6000 N,
  - Back Rest: 3000 N,
  - Maximum load: 180 kg,
  - Maximum Power: 120 VA,
  - Fuse: 1,6 + 6 A,
  - IP44 protection standard,
  - Weight: 90 kg,
  - Base length: 800 mm,
  - Maximum height of the foot rest: 950 - 1000 mm,
  - Sitting height max: 950 mm,
  - Sitting height min: 500 mm,
- **Additional Equipment (extra paid):**
  - Headrest,
  - Paper dispenser,
  - Casters - powder coated metal,
  - Chrome base,
  - IV Pole,
  - Enamelled base,
  - Transparent crystal glass,
  - Base from transparent crystal glass: 19 mm,
- Technical Passport,
- 24-month warranty

