



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Examination Chair TEYCO MED MARTY MT1 (New)

- New,
- Made in Italy,
- Equipped with 2 powerful actuators,
- Electrical Movements:
 - Height adjustment,
 - Back Rest,
- Foot Rest regulated by gas spring,
- Technical Specifications:
 - Power input: 230V - 50 Hz,
 - Power: 24 V,
 - Height: 6000 N,
 - Gas spring: 800 N,
 - Back Rest: 6000 N,
 - Maximum load: 180 kg,
 - Maximum Power: 120 VA,
 - Fuse: 1,6 + 6 A,
 - IP44 protection standard,
 - Waga fotela: 75 kg,
 - Base length: 580 X 740 mm,
 - Sitting height max: 1110 mm,
 - Sitting height min: 650 mm,
- **Additional Equipment (extra paid):**
 - Headrest,
 - IV Pole,
 - 4 gummy wheels with metal frame + breaks,
- Technical Passport,
- 24-month warranty

