



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Stool TEYCO MED WILLY CAR (New)

- New
- Made in Italy
- The perfect blend of innovation, contemporary style and practicality make these doctors stool the best in the market. Aiming to alleviate discomfort and threat of back pain, these stools can make the doctors forget the hours he is putting into work. The practicality and the ergonomic design of the stool also adds to the stool's features
- Vertical movement guaranteed by a gas pump protected by a telescopic chromed cover
- Seat made from polyurethane
- 5 feet basement in chromed steel with footrest ring
- Technical Data:
 - Min height: 490 mm,
 - Max height: 610 mm,
 - Gas pump: 120 mm,
 - Weight: 15 kg,
 - Base diameter: 630 mm,
 - Footrest diameter: 450 mm,
 - Seat diameter: 360 mm,
- Technical Passport
- 24 months warranty
- Financing:
 - Installments
 - Leasing,
 - Loan

