



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medical Display BEACON C22W (New)

- New,
- 2MP Color Medical Display,
- Features:
 - MULTI,
 - BSS,
 - DICOM,
 - SILENCE,
 - VIDEO SWITCH,
 - COLOR,
- Application"
 - Medical Imaging Equipment,
 - Outpatient Service,
 - ICU (Intensive Care Unit),
 - Clinic Department,
 - Community Medical Service,
- Technical Specifications:
 - Panel Size: 21,5",
 - Support Color: 16,7 Billions,
 - Resolution: 1920 x 1080,
 - Pixel Pitch: 0,248 x 0,248 mm,
 - Brightness: 250 cd/m²,
 - Contrast 1000 : 1,
 - Angle of View: 170°,
 - Response Time: 5 ms,
 - Input:
 - DVI-D: 24 pin,
 - DisplayPort: 20 pin,
 - VGA: D-Sub 15 pin,
 - Audio 3,5 mm stereo jack,
 - Speakres: 5W*2,
 - Dimensions: 523 (W) x 376 (H) x 179 (D) mm
 - Weight: 6 kg,
 - Standard Mounting: VESA 100 x 100 mm,
 - Power Supply: AC 230 V / DC 12 V,
 - Max power consumption: 40 W,
 - Power Saving: EPA, VESA,
 - MTBF: > 50 000h,
 - Certificates: EN 60601, CCC, CB, FCC, CE, cTUVus,
- Technical Passport,
- 36-month warranty

