



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Waiting Room Bench (L3TPS1UP) 3 seater (New)**

- New,
- Made in Poland,
- Made from high quality plywood,
- Iron frames painted with powder coating,
- With automatic hinges
- Floor-mounted
- Seats with upholstery
- Dimensions:
  - Width: 168 cm,
  - Depth: 50,5 cm,
  - Height: 80 cm ,
- Available color configurations:
  - Plywood:
    - grey + natural fringes
    - graphite + natural fringes
    - white + natural fringes
    - natural oak
  - Frame:
    - grey,
    - graphite,
    - white
- 24 months warranty

