



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Orthopedic Traction Set ESCHMANN (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in UK,
- Manual,
- Foot support with dual lock,
- Foot support regulated manually,
- Dimensions: 53 x 121 x 72 cm,
- Frame dimensions:
  - 2 x 37,5 cm,
  - 2 x 42 cm,
- Technical Passport,
- 30-days warranty

