



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Orthopedic Traction Set ESCHMANN (Reconditioned)

- Condition: Reconditioned - Very good,
- Made in UK,
- Manual,
- Foot support with dual lock,
- Foot support regulated manually,
- Dimensions: 53 x 121 x 72 cm,
- Frame dimensions:
 - 2 x 37,5 cm,
 - 2 x 42 cm,
- Technical Passport,
- 30-days warranty

