



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



PRAISTECH DEMETER M02 urology chair (New)

New,

Made in Poland,

It has a stable steel construction covered with powder paint,

Additional durability of the chair is ensured by the use of stainless steel

stainless steel is used in the movable parts of the chair,

The backrest consists of three parts, including a removable footrest,

It has a fixed height,

Backrest and seat are adjustable with gas springs, the footrest can be fitted to achieve a lying position,

As standard, the armchair is equipped with knee pads adjustable in all planes, a stainless steel container and a folding platform making it easier to get in and out of the armchair,

There is a possibility to order the upholstery in a colour of your choice, except for "PRESTIGE" colours which are available for an extra charge, Additional options are available (at an extra charge):

Paper hanger,

Caster system,

Hand grips,

Technical parameters:

Length without footrest 139 cm, with footrest 190 cm,

Width 65 cm,

Height 85 cm,

Trendelenburg position 00 - 200,

Anti-Trendelenburg position -150 - 00,

Backrest adjustment angle 00 - 750,

Step height 25 cm,

Weight 70 kg,

Maximum load 150 kg,

It has a valid Technical Passport

Warranty on structure 24 months,

Warranty on hydraulic accessories 12 months,

Possibility of financing (Poland only): Instalments, Leasing, Loan,

If you have any questions, please do not hesitate to contact us!

In case you don't find the product you are interested in, please get in touch with us and we will do our best to find the perfect solution for YOU.

WE SPEAK ENGLISH

HABLAMOS ESPANOL

WIR SPRECHEN DEUTSCH

ON PARLE FRANCAIS

WIJ SPREKEN NEDERLANDS

МЫ ГОВОРИМ ПО РУССКИ

МИ РОЗМОВЛЯЄМО УКРАЇНСЬКОЮ

WE SHIP WORLDWIDE

