



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl


Koordinator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Medicare Standard 2 Section (SM2550) (New)

- New,
- Made in UK,
- Hydraulic Elevation,
- Retract Wheel Base,
- Vertical lifting design with stability throughout,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 194 cm,
 - Width: 63 cm,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from -25° to 75°,
 - Excellent height range from 45 cm to 98 cm,
- **Additional Equipment (extra paid):**
 - Side support rails (Pair) (6002),
 - Side support rails cushions (Pair) (6002P),
 - Breathing hole and plug for prone treatments (6004),
 - Paper Roll Handle (6005),
 - Removable head/knee support (6006),
 - Removable headrest - pillow (6007),
 - Base cover (6053),
 - IV Pole (6054),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD [W CENIE]



TAPICERKI STANDARD - [ZA DOPŁATK.]

