



## Innovation Deluxe 2 Section (MG2675) (New)

- New,
- Made in UK,
- Electric Height Elevation with Electric Backrest and 2-way tilt - linked central locking wheel system,
- Compact frame design,
- 2-way tilt ensures smooth transition when positioning patient for specialist procedures and treatments,
- Basic equipment:
  - Paper roll handle (6005),
  - Removable headrest (6007),
- Anti-bacterial powder coated framework,
- Acrylic capped vac formings to facilitate easy cleaning SEERS ComfortPlus™,
- TiMotion electric linear actuators,
- Table adjustment:
  - Trendelenburg position (from -15° to 15°)
  - Auto-electric CPR functionality,
- Technical Specifications:
  - Length: 200 cm,
  - Width: 70 cm made from 60 mm foam,
  - Safe working load and lifting capacity of 250 kg (528 lbs),
  - Adjustable angle electric backrest from -20° to 85°,
  - Excellent height range from 47 cm to 97 cm,
  - Wheels Tente® Ø125 mm,
- **Additional Equipment (extra paid):**
  - Breathing hole and plug (6004),
  - Side support rails - a pair (6002),
  - Side support cushions - a pair (6002P),
  - Removable head/knee support (6006),
  - Battery back-up (6057),
  - IV Pole (6054),
- Technical Passport,
- 24-month warranty for electric

