



### Medicare Mobile Treatment 3 (SM3650) (New)

- New,
- Made in UK,
- Hydraulic Elevation,
- Adjustable angle of footrest controlled by a gas spring,
- Large Wheel Mobile Base Design for easy movement,
- Compact frame design,
- Basic equipment:
  - Paper roll handle (6005),
  - Side support rails - a pair (6003),
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
  - Length: 191 cm,
  - Width: 63 cm made from foam,
  - Safe working load and lifting capacity of 240 kg (528 lbs),
  - Adjustable angle electric backrest from 0° to 85°,
  - Adjustable angle footrest from - 25° to 75°,
  - Excellent height range from 46 cm to 99 cm,
  - Wheels Tente® Ø125 mm,
- **Additional equipment (extra paid):**
  - Side support cushions - a pair (6003P),
  - Breathing hole and plug (6004),
  - Removable head/knee support (6006),
  - Removable headrest-pillow (6007),
  - Base cover (6053),
  - IV Pole (6054),
- Technical Passport,
- 24-month warranty for electric

