



### **Innovation Standard 3 (MG3670) Section Couch (New)**

- New,
- Made in UK,
- Fully electric height adjustment as well as an electric profiling foot and backrest,
- Individual locking castors or linked central locking wheel system,
- Suitable for a wide range of general examination and minor treatment procedures,
- Acrylic capped vac formings to facilitate easy cleaning,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
  - Length: 200 cm,
  - Width: 65 cm made from 40 mm foam,
  - Safe working load and lifting capacity of 240 kg (528 lbs),
  - Adjustable angle electric backrest from 0° to 80°,
  - Adjustable angle footrest from -35° to 10°,
  - Excellent height range from 45 cm to 95 cm,
  - Wheels - Tente® Ø125 mm,
- **Additional Equipment (extra paid):**
  - Side support rails - a pair (6003),
  - Side support cushions - a pair (6003P),
  - Breathing hole and plug (6004),
  - Paper roll handle (6005),
  - Removable head/knee support (6006),
  - Removable headrest-pillow (6007),
  - IV Pole (6054),
  - Battery back-up (6057),
- Technical Passport,
- 24-month warranty for electric

