



PRAISTON Sp. z o.o.
ul. Górowska 32
64-100 Leszno
Polska

Biuro: +48 65 527 01 67
Serwis: +48 782 844 000
Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
Krzysztof Wybieralski
+48 882 762 006

kw@praiston.pl



Innovation Standard 3 (MG3470) (New)

- New,
- Made in UK,
- Electric Height Elevation with Electric Back and foot rest - Individual locking castors,
- Anti-bacterial powder coated framework,
- Compact frame design,
- Acrylic capped vac formings to facilitate easy cleaning,
- TiMotion electric linear actuators,
- Technical Specifications:
 - Length: 200 cm,
 - Width: 65 cm made from 40 mm foam,
 - Safe working load and lifting capacity of 250 kg (528 lbs),
 - Adjustable angle electric backrest from 0° to 80°,
 - Adjustable angle of footrest fr - 35° to 10°,
 - Excellent height range from 45 cm to 95 cm,
 - Wheels Tente® Ø125 mm,
- **Additional Equipment (extra paid):**
 - Side support rails - pair (6003),
 - Side support cushions - pair (6003P),
 - Breathing hole and plug (6004),
 - Paper Roll handle (6005),
 - Removable head/knee support (6006),
 - Removable headrest - pillow (6007),
 - Battery back-up (6057),
 - IV Pole (6054),
- Technical Passport,
- 24-month warranty for electric

