



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska

Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67

E-mail: biuro@praiston.pl

Koordinator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006

kw@praiston.pl



Innovation Deluxe 3 (MG3475) (New)

- New,
- Made in UK,
- Electric Height Elevation with Electric Back & footrest and 2-way tilt - Individual locking castors,
- 2-way tilt ensures smooth transition when positioning patient for specialist procedures and treatments,
- Basic equipment:
 - Removable headrest (6007)
 - Paper Roll handle (6005),
- Compact frame design,
- Anti-bacterial powder coated framework,
- Acrylic capped vac formings to facilitate easy cleaning SEERS ComfortPlus™,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- Adjustable angle:
 - Trendelenburg position,
 - Anti-Trendelenburg position,
 - Auto-electric CPR functionality and auto-levelling software returns the couch quickly to a fully horizontal position,
- Technical Specifications:
 - Length: 199 cm,
 - Width: 70 cm made from 60 mm foam,
 - Safe working load and lifting capacity of 250 kg (550 lbs),
 - Adjustable angle electric backrest from 0° to 80°,
 - Adjustable angle of footrest from -35° to 10°,
 - Excellent height range from 47 cm to 97 cm,
 - Wheels Tente® Ø125 mm,
- **Additional Equipment (extra paid):**
 - Side support rails - a pair (6003),
 - Side support cushions - a pair (6003P),
 - Breathing hole and plug (6004),
 - Removable head/knee support (6006),
 - Battery back-up (6057),
 - IV Pole (6054),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD (ZA DOPŁATĄ)



