



## Innovation Deluxe Daycare (MG3493) (New)

- New,
- Made in UK,
- Electric Height Elevation with Electric Back & footrest and 2-way tilt - Individual locking castors,
- 2-way tilt ensures smooth transition when positioning patient for specialist procedures and treatments,
- Basic equipment:
  - Removable headrest (6007)
  - Paper Roll handle (6005),
- Compact frame design,
- Anti-bacterial powder coated framework,
- Acrylic capped vac formings to facilitate easy cleaning SEERS ComfortPlus™,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- Adjustable angle:
  - Trendelenburg position,
  - Anti - Trendelenburg position,
  - Auto-electric CPR functionality and auto-levelling software returns the couch quickly to a fully horizontal position,
- Technical Specifications:
  - Length: 199 cm,
  - Width: 70 cm made from 60 mm foam,
  - Safe working load and lifting capacity of 250 kg (550 lbs),
  - Adjustable angle electric backrest from 0° to 80°,
  - Adjustable angle of footrest from -35° to 10°,
  - Excellent height range from 47 cm to 97 cm,
  - Wheels Tente® Ø125 mm,
- **Additional Equipment (extra paid):**
  - Battery back-up (6057),
  - IV Pole (6054),
  - Tray (6241),
- Technical Passport,
- 24-month warranty for electric

