



**PRAISTON Sp. z o.o.**  
 ul. Górowska 32  
 64-100 Leszno  
 Polska

Biuro: +48 65 527 01 67  
 Serwis: +48 782 844 000  
 Fax: +48 65 527 01 67  
 E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
 Krzysztof Wybieralski  
 +48 882 762 006  
[kw@praiston.pl](mailto:kw@praiston.pl)



## **Innovation Deluxe Daycare (MG3493) (New)**

- New,
- Made in UK,
- Electric Height Elevation with Electric Back & footrest and 2-way tilt - Individual locking castors,
- 2-way tilt ensures smooth transition when positioning patient for specialist procedures and treatments,
- Basic equipment:
  - Removable headrest (6007)
  - Paper Roll handle (6005),
- Compact frame design,
- Anti-bacterial powder coated framework,
- Acrylic capped vac formings to facilitate easy cleaning SEERS ComfortPlus™,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- Adjustable angle:
  - Trendelenburg position,
  - Anti - Trendelenburg position,
  - Auto-electric CPR functionality and auto-levelling software returns the couch quickly to a fully horizontal position,
- Technical Specifications:
  - Length: 199 cm,
  - Width: 70 cm made from 60 mm foam,
  - Safe working load and lifting capacity of 250 kg (550 lbs),
  - Adjustable angle electric backrest from 0° to 80°,
  - Adjustable angle of footrest from -35° to 10°,
  - Excellent height range from 47 cm to 97 cm,
  - Wheels Tente® Ø125 mm,
- **Additional Equipment (extra paid):**
  - Battery back-up (6057),
  - IV Pole (6054),
  - Tray (6241),
- Technical Passport,
- 24-month warranty for electric

