



## Examination Table Medicare Multi-Couch - Single Footrest (SM9556) (New)

- New,
- Made in UK,
- Hydraulic Elevation,
- Retractable wheel design with adjustable foot for uneven floors,
- Purpose made precision bearings on all pivot points,
- Clearance beneath frame for mobile hoist,
- Hand switch fitted as standard on electric models,
- Foot sections incorporate anti-trap design,
- TiMotion powerful actuators,
- Technical Data:
  - Length: 180 cm,
  - Width: 65 cm made from foam,
  - Safe working load and lifting capacity (from minimum height) of 240kgs (528 lbs = 37 stone),
  - Adjustable back rest with positive angle of 0 to 80 degrees,
  - Adjustable foot section with angle adjustment of 0 to -85 degrees,
  - Excellent height range from 50 - 105cm,
- **Additional Equipment (extra paid):**
  - Breathing hole and plug for prone treatments (6004),
  - Paper towel holder (6005),
  - Upper limb support (6006),
  - Removable headrest - pillow (6007),
  - Gyneacology stirrups (6044),
  - Gyneacology Leg support (6045),
  - Short base cover (6052),
  - IV Stand (6054),
  - Battery back-up (6057),
  - Single foot switch (6200),
- Technical Passport,
- 24-month warranty for electric, gas spring

