



PRAISTON Sp. z o.o.
 ul. Górowska 32
 64-100 Leszno
 Polska


Biuro: +48 65 527 01 67
 Serwis: +48 782 844 000
 Fax: +48 65 527 01 67
 E-mail: biuro@praiston.pl

Koordynator ds. sprzedaży
 Krzysztof Wybieralski
 +48 882 762 006
kw@praiston.pl



Medicare Orthopaedic-Couches (SM5565) (New)

- New,
- Made in UK,
- Electric,
- Retract Wheel Base,
- Anti-bacterial powder coated framework,
- Vertical lifting design with stability throughout,
- TiMotion electric linear actuators,
- 90 degree dropping extendable foot sections,
- Remote Control,
- Technical Specifications:
 - Length: 180 cm or 185 cm with unfolded footrest,
 - Width: 65 cm made from foam,
 - Safe working load and lifting capacity of 240 kg (528 lbs),
 - Adjustable angle electric backrest from 0° to 80°,
 - Adjustable angle of footrest from 0° to -85°,
 - Excellent height range from 52 cm to 106 cm,
- **Additional Equipment (extra paid):**
 - Orthopaedic arm support (5006),
 - Orthopaedic full leg seat module (5007),
 - Orthopaedic knee support accessory (5008),
 - Orthopaedic leg support attachment (5009),
 - Breathing hole and plug for prone treatments (6004),
 - Removable headrest (6007),
 - Paper Roll handle (6005),
 - Protective Base frame shrouding (6052),
 - Battery Back-up (6057),
 - Operators stool (6103),
 - Deluxe saddle chair (6107),
 - Foot switch (6200),
- Technical Passport,
- 24-month warranty



TAPICERKI STANDARD (W CENIE)



TAPICERKI STANDARD + (ZA DOPŁATĄ)

