



**PRAISTON Sp. z o.o.**  
 ul. Górowska 32  
 64-100 Leszno  
 Polska

Biuro: +48 65 527 01 67  
 Serwis: +48 782 844 000  
 Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
 Krzysztof Wybieralski  
 +48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Medicare Echocardiography Couches (SM3690) (New)


- New,
- Made in UK,
- Electric positioning for height, tilt and backrest operation,
- Adjustable angle of the footrest controlled manually by gas spring,
- Adjustable for Trendelenburg position,
- The large mobile wheel design features central locking and steering facility (LMWD),
- Basic equipment:
  - Removable headrest (6007),
  - Paper roll handle (6005),
- Base frame powder coated with antibacterial layer,
- TiMotion electric linear actuators,
- Remote Control,
- Technical Specifications:
  - Length: 208 cm,
  - Width: 65 cm made from foam,
  - Safe working load and lifting capacity of 240 kg (528 lbs),
  - Adjustable angle electric backrest from 0° to 80°,
  - Adjustable angle footrest from 0° to - 85°,
  - Excellent height range from 50 cm to 103 cm,
- **Additional Equipment (extra paid):**
  - Battery back-up (6057),
  - Patient Lap Tray (6123),
  - Triple foot switch (6209),
- Technical Passport,
- 24-month warranty for electric



Model 6123



50cm 45cm 45cm



Height Range 50cm - 103cm



TAPICERKI STANDARD (W CENIE)

