



**PRAISTON Sp. z o.o.**  
ul. Górowska 32  
64-100 Leszno  
Polska

Biuro: +48 65 527 01 67  
Serwis: +48 782 844 000  
Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordinator ds. sprzedaży  
Krzysztof Wybieralski  
+48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## **Medicare Bariatric 2 Section (SM4542) (New)**

- New,
- Made in UK,
- Bariatric 2 section couch with retractable wheels,
- Remote Control,
- Vertical lifting design with stability throughout,
- Anti-bacterial powder coated framework,
- TiMotion electric linear actuators,
- Technical Specifications:
  - Length: 200 cm,
  - Width: 105 cm,
  - Safe working load and lifting capacity of 325 kg (715 lbs),
  - Adjustable angle electric backrest from 0° to 75°,
  - Excellent height range from 45 cm to 98 cm,
- **Additional Equipment (extra paid):**
  - Breathing hole and plug for prone treatments (6004),
  - Paper Roll handle (6005),
  - Side support rails (Pair) (6074),
  - Side support rails cushions (Pair) (6074P),
  - Memory foam upholstery (60 mm) for additional patient comfort (6076),
  - Double Foot switch (6207),
- Technical Passport,
- 24-month warranty for electric



TAPICERKI STANDARD [W CENIE]



TAPICERKI STANDARD +[ZA DOPŁATKĄ]

