



## Innovation Deluxe Dialysis (MG3690) (New)

- New,
- Made in UK,
- Electric Height Elevation with Electric Back& footrest and 2-way tilt - linked central locking wheel system,
- 2-way tilt ensures smooth transition when positioning patient for specialist procedures and treatments,
- Basic equipment:
  - Removable headrest (6007)
  - Paper Roll handle (6005),
  - Compact frame design,
  - Anti-bacterial powder coated framework,
  - Acrylic capped vac formings to facilitate easy cleaning SEERS ComfortPlus™,
  - Vertical lifting design with stability throughout,
  - Fitted with multi adjustable armrests
  - Fitted with adjustable patient footboard
  - Fitted with pushing handles
  - TiMotion electric linear actuators,
- Adjustable angle:
  - Trendelenburg position,
  - Auto-electric CPR functionality and auto-levelling software returns the couch quickly to a fully horizontal position ,
- Technical Specifications:
  - Length: 199 cm,
  - Width: 70 cm made from 60 mm foam,
  - Safe working load and lifting capacity of 250 kg (550lbs),
  - Adjustable angle electric backrest from 0 ° to 80 °,
  - Adjustable angle of footrest from - 35 ° to 10 °,
  - Excellent height range from 47 cm to 97 cm,
  - Wheels Tente® Ø125mm,
- **Additional Equipment (paid):**
  - Battery back-up (6057),
  - IV Pole (6054),
  - Tray (6241)
- Technical Passport,
- 24 months warranty for electric

