



**PRAISTON Sp. z o.o.**  
 ul. Górowska 32  
 64-100 Leszno  
 Polska

Biuro: +48 65 527 01 67  
 Serwis: +48 782 844 000  
 Fax: +48 65 527 01 67

E-mail: [biuro@praiston.pl](mailto:biuro@praiston.pl)

Koordynator ds. sprzedaży  
 Krzysztof Wybieralski  
 +48 882 762 006

[kw@praiston.pl](mailto:kw@praiston.pl)



## Phlebotomy Couch (SM9556P) (New)

- New
- Made in UK
- Electric with electric backrest and tilt
- Patient handgrip supports - to aid mount/dismount from the couch and provide patient security
- Fully adjustable phlebotomy arms for use on either side of couch
- Fitted with paper roll holder
- Compact frame design with retractable wheel system
- Adjustable foot for stability on uneven surfaces
- Purpose made precision bearings on all pivot points
- The seat section is inclined for patient security
- Clearance beneath frame for mobile hoist
- TiMotion electric linear actuators,
- Technical Data:
  - Length: 180 cm,
  - Standard width 65cm (25")
  - Safe working load and lifting capacity (from minimum height) of 240kgs (528lbs = 37.7 stone)
  - Adjustable angle backrest from -10 to +85 degrees positive
  - Adjustable angle feet plates 0 ° to -85 °,
  - Excellent height range from 50-102cm (19"-39")
- **Additional Equipment (paid):**
  - Removable head cushion (6007)
  - Protective base frame shrouding (6052)
  - High Operators Chair (6101)
  - Deluxe Saddle Seat (6107)
  - Standard Operators Stool (6103),
- Technical Passport
- 24 months warranty for electrics and gas springs



TAPICERKI STANDARD [W CENIE ]



